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Cordiali saluti

RELAB SRL

Raffaello Lobello
Marketing Manager

RELAB SRL
Via C.so Perrone 25 r, 16152 Genova
Telefono 010 667.17.94 Fax 010 612.95.99 email info@relabsrl.it
P.I./C.F 01409780996 C.C.I.A.A Genova R.E.A. n° 413199
Capitale Sociale € 50.000,00 i.v.

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RealCycler®

Validation file

Detection of *Candida auris* DNA
In vitro diagnostic assay

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
Index

1. Introduction.....	3
2. Principle of the test.....	3
3. <i>RealCycler</i> AURI kit components.....	3
4. Validated and compatible purification and amplification systems.....	4
4.1 Purification systems.....	4
4.2 Amplification systems.....	4
5. Specificity.....	5
5.1 Specificity analysis using BLAST algorithm.....	5
5.2 Cross-reactivity evaluation.....	5
6. Analytical sensitivity.....	6
7. Linearity.....	6
8. Validated clinical samples.....	7
9. Stability.....	8
9.1 Stability throughout time assays.....	8
9.2 Freezing/thawing stability assays.....	8
10. Legend symbols.....	9
11. Glossary.....	9

1. Introduction

The present document shows the obtained data in the *RealCycler AURI in vitro* diagnostic reagent validation. It includes the results obtained in specificity studies, analytical sensitivity, linearity and stability studies.

The obtained data prove that the reagent is appropriate for the intended use.

2. Principle of the test

RealCycler AURI is an *in vitro* diagnostic kit of reagents which allows the detection of *Candida auris* DNA in clinical samples. The system includes an internal control of amplification based on the exclusive technology CHIC (*Competitive Heterologous Internal Control*) to prevent false negatives due to reaction inhibition. In Chic-Out *RealCycler* kits, CHIC is out of the AmpliMix to verify the nucleic acid extraction. A Positive Control, which consists of a DNA artificial construction of *Candida auris* is included.

3. *RealCycler AURI* kit components

Component	AURI-M (Monotest)		AURI-U (Universal)		AURI-G (Universal)	
	Vials	Volume	Vials	Volume	Vials	Volume
AmpliMix AURI	24	18,67 µL	1	430 µL	2	430 µL
AURI Positive Control	1	120 µL	1	120 µL	2	120 µL

Component	AURI-UX (Universal)		AURI-GX (Universal)	
	Vials	Volume	Vials	Volume
AmpliMix AURI-UX	1	430 µL	2	430 µL
AURI Positive Control	1	120 µL	2	120 µL
CHIC	1	350 µL	2	350 µL

4. Validated and compatible purification and amplification systems

4.1 Purification systems

Purification systems	
Validated for AURI / -U / -G	Compatible for AURI / -U / -G
<p>QIAamp DNA Blood Mini kit (references 51104, 51106). Qiagen.</p>	<p>Arrow/LIAISON[®] Ixt. DiaSorin BioRobot EZ1. Qiagen[®] QIAcube[®]. Qiagen[®] MagCore[®] Automated Nucleic Acid Extractor. RBCBioscience. Maxwell[®] 16 Cell LEV DNA Purification kit. Promega Corporation. NucliSENS[®] eas, MAG[®]. bioMérieux.</p>
Validated for AURI / -U / -G	Compatible for AURI / -U / -G
<p>QIAamp Viral RNA Mini Kit (references 52904, 52906). Qiagen. MagCore[®] Automated Nucleic Acid Extractor (references MVN400-03, MVN400-04). RBCBioscience. Maxwell[®] 16 Vial Total Purification kit (reference AS1150). Promega Corporation.</p>	<p>Arrow/LIAISON[®] Ixt. DiaSorin BioRobot Ez1. Qiagen[®] QIAamp DNA Blood Mini Kit. Qiagen[®] QIAcube[®]. Qiagen[®] Maxwell[®] 16 Cell LEV DNA Purification kit. Promega Corporation. NucliSENS[®] easyMAG[®]. bioMérieux.</p>

4.2 Amplification systems

Validated Real-time PCR Instruments
<p>SmartCycler[®]. Cepheid[®] CE896[®]. Bio-Rad[®] ABI 7500. Thermo Fisher Scientific Rotor Gene[®] Q. Qiagen[®] LightCycler[®] 480. Roche Diagnostics</p>

5. Specificity

5.1 Specificity analysis using BLAST algorithm

The analysis using BLAST algorithm (www.ncbi.nlm.nih.gov/blast) shows that selected primers and probes detect *Candida auris* Targets used in pathogens detection:

Pathogen	Target
<i>Candida auris</i>	Internal Transcribed Spacer Region

5.2 Cross-reactivity evaluation

Pathogen	RealCycler AURI
	<i>Candida auris</i>
<i>Candida albicans</i>	Not detected
<i>Candida glabrata</i>	Not detected
<i>Candida krusei</i>	Not detected
<i>Candida parapsilosis</i>	Not detected
<i>Candida tropicalis</i>	Not detected
<i>Neisseria meningitidis</i>	Not detected
<i>Staphylococcus aureus</i>	Not detected
Human genomic DNA	Not detected

6. Analytical sensitivity

The analytical sensitivity limit of *RealCycler* AURI / -U / -G / -UX / -GX kit has been determined by limit dilution using a CFX96™ (Bio-Rad) instrument. This value is 1 copy/μL for *Candida auris*.

The analytical sensitivity has been tested on formats which include CHIC in the AmpliMix. In Chic-Out formats, the sensitivity will depend on the extraction system efficacy.


7. Linearity

It has been tested on formats which include CHIC in the AmpliMix. In Chic-Out formats, the linearity would depend on the extraction system efficacy.

Ct value (Cycle Threshold) obtained in each pathogen dilution.

Pathogen	Concentration (copies/μL)								
	10 ⁸	10 ⁷	10 ⁶	10 ⁵	10 ⁴	10 ³	10 ²	10 ¹	1
<i>Candida auris</i>	8,31	12,09	15,29	18,78	22,15	25,40	29,09	32,21	36,50

Linearity in *Candida auris* detection

In the range 1 copy/μL – 10⁸ copies/μL, the technique maintains the linearity (r²=0,9993).

8. Validated clinical samples

RealCycler AURI has been validated with the following clinical samples:

Pathogen	Sample
<i>Candida auris</i>	Blood culture

9. Stability

9.1 Stability throughout time assays

The objective of this study is to evaluate the reagent stability throughout time, and to verify that it preserves the specifications up to the expiration date, under the recommended conservation conditions. This study compares the Ct value obtained with the evaluated reagent, with the Ct value obtained by a previous validated reference reagent.

Results obtained in stability throughout time assays:

Candida auris

Time	Expected Ct	Obtained Ct
Month 1	30,65	30,73
Month 3	30,65	30,90
Month 6	30,65	30,61
Month 9	30,65	30,49
Month 12	30,65	31,14
Month 18	30,65	31,26

CHIC

Time	Expected Ct	Obtained Ct
Month 1	33,98	33,98
Month 4	33,98	34,03
Month 6	33,98	33,68
Month 9	33,98	34,68
Month 12	33,98	33,45
Month 18	33,98	33,54

Expected Ct: Ct value obtained with a reference AmpliMix lot.

Obtained Ct: Ct value obtained with the AmpliMix subjected to stability studies.


Results show that Ct values obtained in the *Candida auris* and CHIC study do not show significant variations during the study period.

9.2 Freezing/thawing stability assays

The objective of this study is to evaluate the reagent stability after freezing/thawing cycles. This study compares the Ct value obtained with the evaluated reagent during the freezing/thawing cycles assay.

The experimental assays performed have shown that *RealCycler* reagents are stable for at least 15 cycles of freezing / thawing.


10. Legend symbols

Symbol	Meaning
	Manufacturer
	Catalog number
	Registered as <i>in vitro</i> diagnostic reagent in the European community
	Reaction number

11. Glossary

Ct: Cycle Threshold. The cycle of the PCR protocol in which appears significant fluorescence is proportional to the DNA/RNA quantity present in the sample.

Guidance for Detection of Colonization of *Candida auris*

Candida auris is an emerging multidrug-resistant yeast that can cause invasive infections with high mortality and has been transmitted in healthcare settings. Patients may be asymptotically colonized with *C. auris*. Patients with *C. auris* colonization can spread this yeast to other patients, and colonized patients can develop invasive as well as superficial infections. Identifying persons colonized with *C. auris* is a key step in containing the spread of *C. auris*. The instructions below will provide guidance for processing swabs to assess for *C. auris* colonization.

For this testing, the currently recommended sites are axilla, groin, and sometimes nares. For more information on the screening of patients, including a procedure for the collection of patient swabs, fact sheets for patients, and a sample script for talking with patients about screening, read CDC's information on screening for laboratory staff and health professionals.

For additional information regarding *C. auris*, including the latest identification, treatment, and infection control recommendations, visit CDC's *C. auris* home page.

All confirmed *C. auris* specimens should be reported to state and local health departments and CDC at candidaauris@cdc.gov. Laboratories are encouraged to submit *C. auris* isolates to the Antibiotic Resistance (AR) Lab Network.

Disclaimer

The techniques for isolation presented here have not been completely validated for diagnostic purposes. Any laboratory wishing to employ these procedures should complete a validation that conforms with their regulatory oversight. None of these tests have been cleared or approved by the FDA. Mention of a specific commercial product does not qualify as an endorsement of the product by the Centers for Disease Control and Prevention.

Real-Time PCR

For those institutions with the resources to do so, real-time PCR is a fast and accurate method for detecting *C. auris*. Institutions may use this protocol [\[PDF – 11 pages\]](#) for screening patients using real-time PCR. There are also other publications that provide protocols for PCR-based identification of *C. auris*.

Alternatives to PCR

For institutions without the resources to perform real-time PCR, a number of options are available for the screening of patients using culture-based methods. This procedure can be optimized to fit with current laboratory capacity. The following outlines the steps to isolate *C. auris* from patient swabs.

Suggested equipment/material for culture-based methods

- Culture collection and transport system that can be used to swab a patient. Examples include:
 - Nylon-flocked swab (BD Eswab collection and transport system; Becton Dickinson and Company, Sparks, MD)
 - Rayon tip swabs (Fisherfinest Amies Charcoal bacteriology culture collection and transport system; Fisher Healthcare, Ontario, Canada)
- CHROMagar *Candida* Chromogenic agar
 - There are multiple commercial companies that sell CHROMagar *Candida* media as well as at least one company that sells powder for making the media in-house
- 10 µl transfer loops
- Vortex
- Stationary Incubator

Considerations

Communities of *Candida*

Swabs from patients will contain a complex community of microorganisms including species of *Candida* and other yeasts. It is necessary to isolate *C. auris* from this community to obtain accurate species identification.

***Note:** All the following culture methods include an incubation step to help grow and isolate *C. auris* from this community. Incubation should optimally take place at 40°C. *Candida auris* can grow at temperatures up to 42°C, but 37°C–40°C is optimal to account for temperature fluctuations in the incubator. If a 40°C incubator is not available, 37°C will work.

Strategies for screening non-sterile body sites

There are many methods for processing patient specimens to look for *C. auris* colonization. This list is not meant to be comprehensive but to give laboratories several options that may be incorporated into their current laboratory workflow.

Salt/Dulcitol enrichment broth

This broth is commercially available from S2 media (Spokane Valley, WA) and Thermo Oxoid (Ontario, Canada), but the recipe has been published here [7]. This broth contains 10% NaCl, which inhibits the growth of most, but not all, yeast species. *Candida glabrata* is one of the species that can sometimes grow at high salt, and therefore dulcitol is used as a carbon source rather than glucose. *C. glabrata* cannot assimilate or ferment dulcitol. Although CDC has primarily used this media in broth form, it could also be used as a salt/dulcitol agar plate.

This media is selective but not differential. A differential media such as CHROMagar will still have to be used to differentiate between *C. auris* and other species that grow on the plates.

Plating directly to salt/dulcitol agar

There are a couple methods for direct plating. The first is to remove the swab and to roll it over the agar plate, covering all surfaces. The second option is to vortex and disperse the pathogen from the swabs before removing liquid from the transport, transfer it to the plate, and then use a spreader to cover all surfaces. The advantage of using liquid from the transport media is that a larger volume can be plated, increasing your chance of recovering *C. auris* from a patient who is not heavily colonized. The disadvantage is that it may be difficult to pick a single colony from a patient who is highly colonized or who may be colonized by two or more species of yeast.

Again, this media is selective but not differential. A differential media such as CHROMagar will still have to be used to differentiate between *C. auris* and other species that grow on the plates.

Direct plating may be faster and easier for many laboratories than the salt/dulcitol enrichment broth procedure. However, some laboratories have found direct plating to be less sensitive than the enrichment broth method.

Agar media CHROMagar *Candida*

There are no conventional agar plates that allow for the differentiation of *C. auris*.

Chromagar *Candida* allows the differentiation between *C. albicans*, *C. tropicalis*, *C. krusei*, and everything else. On CHROMagar, *C. auris* can be differentiated from the most common species, including *C. glabrata*. However, *C. auris* does not have a single distinct color on CHROMagar; *C. auris* can be pink, red, cream, and even purple, as well as a combination of these colors on the same plate (see images). While the manufacturer's recommended incubation temperature for CHROMagar is 37°C, we have found that incubation at 40°C will inhibit the growth of some other organisms.

This means that all suspect colonies must be identified to species. MALDI-TOF (Bruker Biotyper or bioMérieux VITEK MS), DNA sequencing of marker genes D1/D2 or the ITS region, or whole genome sequencing provide the most accurate species identification.

If using phenotypic methods, an algorithm [PDF - 11 pages] can be used to identify *C. auris*.

Safety considerations


If possible, work with *C. auris* should be performed in a biological safety cabinet or a glove box to ensure safety and sterility.
<https://www.cdc.gov/fungal/candida-auris/c-auris-guidance.html>

9/3/2021

If possible, work with *C. auris* should be performed in a biological safety cabinet or a glove box to ensure safety and sterility.

C. auris can survive for weeks on surfaces, has reduced susceptibility to quaternary ammonia disinfectants, and can colonize skin. Therefore, strict BSL2 laboratory safety precautions must be followed when working with this organism. Specifically, it is recommended that cultures are processed within BSL2 biosafety cabinet, gloves and lab coats should be worn, and strong hand hygiene should be employed. Since products with *C. albicans* or fungicidal claims may not be effective against *C. auris*, 10% bleach or products with EPA approval for *C. auris* should be used for cleaning the work area. A secondary container should be used when transporting *C. auris*, and this container should be disinfected after use.

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Trattamenti e il copyright sono: (0102-196-2602)